→ BPL_TEST2_Batch script with FMPy ver 0.3.15

The key library FMPy ver 0.3.15 is installed.

After the installation a small application BPL_TEST2_Batch is loaded and run. You can continue with this example if you like.

```
!lsb release -a # Actual VM Ubuntu version used by Google
    No LSB modules are available.
    Distributor ID: Ubuntu
    Description: Ubuntu 20.04.5 LTS
                   20.04
    Release:
    Codename:
                   focal
%env PYTHONPATH=
    env: PYTHONPATH=
!wget https://repo.anaconda.com/miniconda/Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
!chmod +x Miniconda3-py39_23.1.0-1-Linux-x86_64.sh
!bash ./Miniconda3-py39_23.1.0-1-Linux-x86_64.sh -b -f -p /usr/local
sys.path.append('/usr/local/lib/python3.9/site-packages/')
    --2023-06-01 06:42:02-- https://repo.anaconda.com/miniconda/Miniconda3-py39_23.1.0-1-Linux-x
    Resolving repo.anaconda.com (repo.anaconda.com)... 104.16.130.3, 104.16.131.3, 2606:4700::6810
    Connecting to repo.anaconda.com (repo.anaconda.com) | 104.16.130.3 | :443... connected.
    HTTP request sent, awaiting response... 200 \mbox{OK}
    Length: 69888122 (67M) [application/x-sh]
    Saving to: 'Miniconda3-py39 23.1.0-1-Linux-x86 64.sh'
    Miniconda3-py39 23. 100%[===========] 66.65M
                                                             149MB/s
    2023-06-01 06:42:02 (149 MB/s) - 'Miniconda3-py39 23.1.0-1-Linux-x86 64.sh' saved [69888122/69
    PREFIX=/usr/local
    Unpacking payload ...
    Installing base environment...
    Downloading and Extracting Packages
    Downloading and Extracting Packages
    Preparing transaction: done
    Executing transaction: done
    installation finished.
!conda update -n base -c defaults conda --yes
```

!conda --version
!python --version

conda 23.5.0 Python 3.9.16

!conda install -c conda-forge fmpy --yes # Install the key package

Preparing transaction: done Verifying transaction: done Executing transaction: done

!conda install matplotlib --yes

```
Preparing transaction: done
Verifying transaction: done
Executing transaction: done

#!conda install scipy --yes

#!conda install xlrd --yes

#!conda install openpyxl --yes
```

→ BPL_TEST2_Batch setup

Now specific installation and the run simulations. Start with connecting to Github. Then upload the two files:


- FMU BPL_TEST2_Batch_linux_om_me.fmu
- Setup-file BPL_TEST2_Batch_fmpy_explore.py

```
Linux - run FMU pre-compiled OpenModelica 1.21.0
    Model for bioreactor has been setup. Key commands:
     - par()
               - change of parameters and initial values
     - init()
                  - change initial values only

    simulate and plot


     - simu()
     - newplot() - make a new plot
                  - show plot from previous simulation
     - show()
                   - display parameters and initial values from the last simulation
     - disp()
     - describe() - describe culture, broth, parameters, variables with values/units
    Note that both disp() and describe() takes values from the last simulation
    Brief information about a command by help(), eg help(simu)
    Key system information is listed with the command system_info()
%matplotlib inline
plt.rcParams['figure.figsize'] = [25/2.54, 20/2.54]
import warnings
warnings.filterwarnings("ignore")
```

▼ BPL_TEST2_Batch - demo




Simulation were initial value of biomass VX_0 is varied
newplot(plotType='TimeSeries')
for value in [1.0, 0.7, 1.3]: init(VX 0=value); simu(5)

Restore default value of VX_0
init(VX_0=1.0)



```
# Simulation were initial value of substrate VS_0 is varied
newplot(plotType='TimeSeries')
for value in [10, 7, 13]: init(VS_0=value); simu(5)


# Restore default value of VS_0
init(VS_0=10)
```



```
# Simulation where metabolism is changed after 3 hours
newplot(plotType='TimeSeries')
simu(5)

simu(3)
par(Y=0.4, qSmax=1.0/(0.4/0.5)); simu(2, 'cont')

# Restore default value of Y and qSmax
par(Y=0.5, qSmax=1.0)
```



√ 0s completed at 08:48